

The Porter Midas™ Digital Sedation Machine

The Superior Choice in Modern Sedation Technology

Feature Comparison

Feature	Midas	Accutron Digital	Matrix Digital
Touchscreen Interface	YES	NO	NO
Accuracy (under 5LPM & Above 5LPM)	+/- 0.3% & 0.05	+/- 0.5% & 0.05	+/- 0.5% & 0.05
Bluetooth & Data Log	YES	NO	NO
Nitrous Oxide Lock	YES	NO	NO
Gas Usage Tracking	YES	NO	NO



[Click to watch Video](#)



Only with Midas Digital!

- **Full Touchscreen Control:** Intuitive digital interface for precise, fast adjustments—industry-leading usability. Works even with gloves on!
- **Electronic Vacuum Scavenging:** Real-time monitoring & automated eAVS control for consistent safety.
- **Session Log & Bluetooth Connection:** Digital session logging for compliance as well as bluetooth connection to iPad app.
- **Nitrous Lock:** Pin activated Nitrous oxide lock to reduce misuse and incorrect gas delivery.
- **Comprehensive Safety Suite:** Includes fail-safes, alarms, and automatic emergency cut-offs for peace of mind.
- **Gas Usage Tracking:** Tracking of gas use including total litres dispensed and percentage remaining.
- **Session Timer & Clock Display:** On-screen timer for session duration and clock display.

Other Key Features

- Adjustable Gas Concentration
- Remote Cabinet / Trolley Mount
- Safety Alarm System
- DISS Gas Fittings
- Oxygen Therapy and Emergency Oxygen Flush setting