

YUWELL BREATHWEAR MASK SERIES



Fit



Comfort



Convenience



Light & Soft




Rebound Silicone Cushion

The mask is made of rebound silicone cushion, which is airtight, light and soft. It is designed to reduce facial pressure to make it fit and comfortable.



Adaptive Silicone Cushion

The silicone cushion has an adaptive design at the nose bridge, which can be automatically adjusted according to different face shapes, so as to better relieve the pressure on the high nose and improve user experience.



Quiet Vent-Hole

Honeycomb design to reduce airflow erosion intensity. Gentle and quiet.



Quick-Release Connector

Using a quick release connector to easily separate or connect the mask and the device.





Soft and Light Headgear

Selected soft skin-friendly fabrics, light and stable, help patients adapt to the treatment comfortably.



| TECHNICAL SPECIFICATIONS | | | | |
|--|---|---|---|--|
| Picture |  |  |  |  |
| Model | YF-01 | YN-02 | YF-03 | YN-01 |
| Certification | CE&FDA | CE&FDA | CE&FDA | CE&FDA |
| Type | Full Face Mask | Nasal Mask | Full Face Mask | Nasal Mask |
| Therapy Pressure Range | 4-30 cmH ₂ O | 4-30 cmH ₂ O | 4-30 cmH ₂ O | 4-30 cmH ₂ O |
| Size | S/M/L/XL | S/M/L | S/M/L | S/M/L |
| Tube Connection | φ22mm | φ22mm | φ22mm | φ22mm |
| Flow Resistance(At 50 L/min) | ≤0.5cmH ₂ O | ≤0.5cmH ₂ O | ≤0.5cmH ₂ O | ≤0.5cmH ₂ O |
| Flow Resistance(At 100 L/min) | ≤1.5cmH ₂ O | ≤1.5cmH ₂ O | ≤1.5cmH ₂ O | ≤1.5cmH ₂ O |
| Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min) | ≤2.0cmH ₂ O/(L/S) | — | ≤2.0cmH ₂ O/(L/S) | — |
| Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min) | ≤2.0cmH ₂ O/(L/S) | — | ≤2.0cmH ₂ O/(L/S) | — |
| Anti Asphyxia Valve Open Pressure | ≤2.0cmH ₂ O/(L/S) | — | ≤2.0cmH ₂ O/(L/S) | — |
| Anti Asphyxia Valve Closed Pressure | ≤2.0cmH ₂ O/(L/S) | — | ≤2.0cmH ₂ O/(L/S) | — |
| G.W./N.W. | 147/106g | 90/60g | 139/103g | 101/63g |
| Packing | 10PCS/CTN 35×28×28cm | 10PCS/CTN 35×24×24cm | 10PCS/CTN 35×28×28cm | 10PCS/CTN 35×24×24cm |

| TECHNICAL SPECIFICATIONS | | | |
|--|---|---|---|
| Picture |  |  |  |
| Model | YF-02 | YN-03 | YP-01 |
| Certification | CE&FDA | CE&FDA | CE&FDA |
| Type | Full Face Mask | Nasal Mask | Pillow Mask |
| Therapy Pressure Range | 4-30 cmH ₂ O | 4-30 cmH ₂ O | 4-20 cmH ₂ O |
| Size | S/M/L/XL | S/M/L | S/M/L |
| Tube Connection | φ22mm | φ22mm | φ22mm |
| Flow Resistance(At 50 L/min) | ≤0.5cmH ₂ O | ≤0.5cmH ₂ O | ≤0.5cmH ₂ O |
| Flow Resistance(At 100 L/min) | ≤1.5cmH ₂ O | ≤1.5cmH ₂ O | ≤1.5cmH ₂ O |
| Expiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min) | ≤2.0cmH ₂ O/(L/S) | — | — |
| Inspiratory Resistance with Anti Asphyxia Valve Open (At 50 L/min) | ≤2.0cmH ₂ O/(L/S) | — | — |
| Anti Asphyxia Valve Open Pressure | ≤2.0cmH ₂ O/(L/S) | — | — |
| Anti Asphyxia Valve Closed Pressure | ≤2.0cmH ₂ O/(L/S) | — | — |
| G.W./N.W. | 148/107g | 129/99g | 268/81g |
| Packing | 10PCS/CTN 35×28×28cm | 10PCS/CTN 35×28×28cm | 10PCS/CTN 34×22×24.5cm |