



**MEGA
MEDICAL**

www.megamedical.com.au

**ANALGESIC AND ANAESTHESIA +
PENTHROX**

Analgesic and Anaesthesia	1
Sedation Equipment	1
Sedation Nasal Hoods and Consumables	5
Analgesic and Anaesthesia Mounts.....	8
Penthrox	9

ANALGESIC AND ANAESTHESIA

SEDATION EQUIPMENT



PORTER MXR SEDATION FLOWMETER PACKAGE

Brand: Porter

SKU: C3400A

Can also be purchased with the Low-profile Silhouette Mask system

Description

Porter's anaesthesia equipment is designed with Australia's clinicians and medical teams in mind, giving them the capability required to achieve excellent levels of care for patients. This package from Mega Medical includes: MXR Series Nitrous Oxide Sedation Flowmeter and Porter's autoclavable patient breathing circuit. This unit is fully compliant with all Australian standards including AS2896 and AS3840 and also complies to the ADA Policy 6.17 for Conscious Sedation.

The MXR-3000 has been the most trusted and recognised flowmeter on the market, featuring Porter precision glass indicator tubes and ergonomic single motion controls. The MXR 3000 is designed for those who would prefer an at-your-side vertical flowmeter.

Another superior feature of the Porter MXR Series is the control it provides to clinicians during patient care. Clinicians can set the O₂ and N₂O mix ratio by regulating the N₂O percentage. This enables effective sedation and anaesthesia for all patients, including paediatric care cases. The controls of the device are also more responsive than other units on the market, enabling swift dose adjustments as and when required.

Mega Medical is the licensed distributor for these products on the Australian market. ARTG Entry Number: 375145

Features

Porter MXR 3000

- On/off switch and oxygen flush button
- The Porter MXR can be equipped with either a standard breathing circuit with In-Line vacuum control or a breathing circuit with the patented Automatic Vacuum Switch (AVS) control.
- A clearly defined scale provides better ease of use than a digital display.
- No power is required to run analog flowmeter system.
- Easier to service and repair compared to digital models.

Standard Contents

- The Porter MXR Flowmeter
- 3 Litre Silicone Reservoir Bag
- 1 x Scavenging circuit with in-line Vacuum Control
- 1 x Autoclavable mask and 3 x liners (Adult)
- Nitrous and Oxygen Elbow adapters to comply with Australian Standards

*Silhouette Mask (Optional Accessory) - More info pg9

- Low profile design with tapered fresh gas outlet to fit perfectly over the patient's airways.
- Improved infection control in line with hospital / medical standards.
- Clear mask allows a better observation of patients breathing.
- Lightweight, compact and easy to store.
- Single use.



Scavenging System

- Unique design, comfortable, quiet and safe.
- Features an outer hood and the option for either a flexible autoclavable or single-use-disposable inner liner.
- Includes 3 litre reservoir bag, fresh gas tubing, scavenger breathing coaxial tubing, adult or paediatric autoclavable hood (with 3 inner liners), vacuum line hose and options to include the Porter Automatic Vacuum Switch or In-Line Vacuum Control.
- The inner liner creates a complete, yet comfortable seal around the patient's nose – thus assuring the Porter mask's superior scavenging efficiency.

Related Products



Porter 5 Inch
Cabinet Mount
SKU: 2035-1



Porter MXR Tele-
scoping Wall Mount
SKU: 2020



Porter Autoclavable
Hood Kit – Adult
SKU: 5054A



Porter Autoclavable
Hood Kit – Paediatric
SKU: 5054B



Porter Sedation
Breathing Circuit –
Adult
SKU: 5155-2



Porter Sedation
Breathing Circuit –
Paediatric
SKU: 5155-1



3 Litre Sedation
Rebreathing Bag
SKU: AMBBT1



1 Litre Sedation
Rebreathing Bag
SKU: PA-50259



Corrugated Fresh Gas
Tubing – Nitrous Oxide
Flowmeter
SKU: AMBCT1



Tri Cylinder Trolley
SKU: AMBST01



Porter Mount for Tri
Cylinder Trolley
SKU: POR4527



PORTABLE SEDATION SYSTEM PACKAGE – PORTER MXR + AUTOCLAVABLE BREATHING CIRCUIT

Brand: Mega Medical

SKU: meg-C3400A1

Can also be purchased with the Low-profile Silhouette Mask system

Description

Porter's anaesthesia equipment is designed with Australia's clinicians and medical teams in mind, giving them the capability required to achieve excellent levels of care for patients. This package from Mega Medical includes: MXR Series Nitrous Oxide Sedation Flowmeter and porter's autoclavable patient breathing circuit. This unit is fully compliant with all Australian standards including AS2896 and AS3840 and also complies to the ADA Policy 6.17 for Conscious Sedation.

The MXR-3000 has been the most trusted and recognised flow meter on the market for a number of years. Featuring Porter precision glass indicator tubes and ergonomic single motion controls it is the most user-friendly unit available. The MXR 3000 is designed for those who would prefer an at-your-side vertical flow meter.

Another superior feature of the Porter MXR Series is the control it provides to clinicians during patient care. Clinicians can set the O₂ and N₂O mix ratio by regulating the N₂O percentage. This enables effective sedation and anaesthesia for all patients, including paediatric care cases. The controls of the device are also more responsive than other units on the market, enabling swift dose adjustments as and when required.

Standard Contents

- Porter MXR Flowmeter
- Adult Autoclavable Breathing Circuit
- Tri Cylinder Trolley
- Medireg Oxygen Regulator
- Medireg Nitrous Oxide Regulator
- Oxygen Hose Assembly
- Nitrous Oxide Hose Assembly
- Flowmeter Trolley Mount

Features

Porter MXR 3000

- On/off switch and oxygen flush button
- The Porter MXR can be equipped with either a standard breathing circuit with In-Line vacuum control or a breathing circuit with the patented Automatic Vacuum Switch (AVS) control
- A clearly defined scale provides better ease of use than a digital display
- No power is required to run analog flowmeter system
- Easier to service and repair compared to digital models

*Silhouette Mask (Optional Accessory) - More info pg9

- Low profile design with tapered fresh gas outlet to fit perfectly over the patient's airways.
- Improved infection control in line with hospital / medical standards.
- Clear mask allows a better observation of patients breathing.
- Lightweight, compact and easy to store.
- Single use.





ENTONOX DEMAND VALVE + 1M HOSE AND MASK

Brand: Mega Medical

SKU: MEGA-EQ

Description

The Mega Entonox (Equanox) administration system is a rugged demand inhalation unit for the administration of 50% nitrous oxide / 50% oxygen analgesic gas mixture (Entonox or Equanox) and was specifically designed for use in emergency and obstetric analgesia.

The system consists of a single-stage pressure-reducing regulator, a reinforced supply hose, a demand valve and a universal face mask and bleed pin.

Standard Contents

- Demand Valve – Entonox (DVR)
- 1m Hose Assembly
- LR5134 Mask

Specifications

- TGA Registered – #225911
- Input Pressure: 50 PSI
- Operating Temperature: -18 to 50°C
- Demand Breathing Flow Rate: 0-140 L/min
- Storage Temperature: -40 to 60°C
- Patient Valve Dead Space: 8 ml
- Relative Humidity for Storage and Operating Use: 15 to 95%
- Dimensions: Diameter (mm) 57mm x 69mm (approx.)
- Weight: 200g (approx.)

Related Products




Entonox Hose Assembly – Demand fitting – 1m
SKU: CIGMA167



Analgesic Inhalation System
SKU: MEGA-EQ550010



Disposable Mouth-piece 22mm
SKU: PH56001



Disposable Mouth-piece - 22mm - BOX50
SKU: PH56001-BOX50



Entonox Demand Valve
SKU: AMB-EQ



ADULT FRESH MINT – CLEARVIEW SEDATION HOODS (BOX 12)

Brand: Accutron

SKU: GDAC33035-16FM

Description

The Accutron ClearView single-use nasal hood is used for the delivery of nitrous oxide conscious sedation to a patient. The innovative double hood design from Accutron further reduces ambient nitrous oxide in the operating room, while the clear outer hood provides a clear visual to the clinician when a patient breathes through their nose (condensation appears in hood).

Features

- Hood-in-hood design draws excess gas into vacuum reducing ambient nitrous oxide in operating room
- Clear outer hood shows patient breathing through nose respiration (condensation appears in hood)
- Soft inner hood creates desirable nasal/facial seal
- Low profile expands clinical field of view
- Appealing colours and scents engage and relax patients
- Single-use minimises cross-contamination; saves staff time
- Latex free material protects patients and staff who have latex allergies

Variations



Adult French Vanilla
– ClearView Sedation Hoods (box 12)
SKU: GDAC33035-14FV



Adult Unscented Grey
– ClearView Sedation Hoods (box 12)
SKU: GDAC33035-9UG



Adult Large Unscented Grey
– ClearView Sedation Hoods (box 12)
SKU: GDAC33034-9M



Paediatric Outlaw Orange – ClearView Sedation Hoods (box 12)
SKU: GDAC33037-1000



Paediatric Birthday Bubblegum – ClearView Sedation Hoods (box 12)
SKU: GDAC33037-12BB



Paediatric Strawberry – ClearView Sedation Hoods (box 12)
SKU: GDAC33037-16SS



Paediatric Grape – ClearView Sedation Hoods (box 12)
SKU: GDAC33037-17G



Paediatric Unscented Grey – ClearView Sedation Hoods (box 12)
SKU: GDAC33037-9UG



SILHOUETTE VARIETY PACK – 4 SIZES X 3EA

Brand: Parker Porter

SKU: SIL-VAR-4x3



Description

The Silhouette nasal mask and breathing circuit are designed with lightweight tubing and conforming mask used to deliver nitrous oxide and oxygen to a conscious, spontaneously breathing patient. The Silhouette nasal mask also provides scavenging capability when connected to a vacuum source to remove excess nitrous oxide from the patient's exhalation.

Features

- Low profile
- Transparent sedation hoods to monitor patients breathing more effectively
- High accessibility to the patient oral cavity
- Lower infection rate
- Improved patient experience and comfort
- Effective scavenging system

Specifications

- 4 sizes, 3 masks ea (Paediatric, Small, Medium, Large)
- Latex free
- Low profile
- Clear sedation masks, available in 4 sizes and colour coded
- Perfect seal and high accessibility to the patient oral cavity
- Reduced infection rate during procedures

Related Products



PORTER CONVERSION KIT FOR SILHOUETTE MASK SYSTEM

Brand: Parker Porter

SKU: SIL-Conn-Kit

Description

The Silhouette conversion kit will allow you to convert any analog sedation unit to work with the Silhouette Sedation Mask System. The Porter Silhouette Mask System provides revolutionary nitrous oxide delivery to the patient in a low profile completely disposable setup while providing effective scavenging ability at the same time.

Specifications

- 200cm white flexible fresh gas and scavenging tube
- Fresh gas tube – 6mm ID, Scavenging tube – 10mm ID



SEDATION NASAL HOODS AND CONSUMABLES



Porter Adult Disposable Liners
– Pack 12 – Unscented
SKU: 5053 AD12



Porter Paediatric Disposable Liners
- 12 Pack - Unscented
SKU: 5053 PD12



Paediatric Disposable Liners - Pack
12 - Citrus
SKU: 5053 PD12C



Autoclavable Nasal Hood – Large
SKU: GDAC33018LG



Scavenging Tubing – 1m
SKU: 6789912

Variations



Autoclavable Nasal Hood – Medium
SKU: GDAC33019Med



Autoclavable Nasal Hood – Small
SKU: GDAC33020Sm

ANALGESIC AND ANAESTHESIA MOUNTS



Wall Mounting Stand – Sedation Mixer
SKU: AMBST02



Bi-Fold Wall Arm
SKU: GDAC32001



Single Articulating Folding Arm
SKU: GDAC32002



Telescoping Wall Mount
SKU: GDAC32003



Porter 5 Inch Cabinet Mount
SKU: 2035-1



Porter MXR Telescoping Wall Mount
SKU: 2020

Penthrox



PENTHROX SINGLE COMBO INHALER

Brand: MDI

SKU: MDI-ME-MS245

Description

Penthrox (methoxyflurane) is a non-narcotic inhalable analgesic that is self-administered by a conscious haemodynamically stable patient (under supervision of personnel trained in its use) via the hand-held Penthrox Inhaler to provide pain relief.

Features

- It is clear
- Almost colourless liquid with a characteristic fruity smell that becomes a vapour or gas once it has been poured into the Penthrox Inhaler



Standard Contents

- 1x 3ml Penthrox (methoxyflurane)
- 1x Penthrox Inhaler

Variations



Penthrox Single Combo Inhaler – Pack 10

SKU: MDI-ME-MS245-10

PENTHROX SINGLE COMBO WITH CHAMBER

Brand: MDI

SKU: MDI-ME-MS260

Description

Penthrox (methoxyflurane) is a non-narcotic inhalable analgesic that is self-administered by a conscious, haemodynamically stable patient (under supervision of personnel trained in its use) via the hand-held Penthrox Inhaler to provide pain relief.

Standard Contents

- 1x 3ml Penthrox (methoxyflurane)
- 1x Penthrox Inhaler
- 1x AC Chamber

Features

- It is clear
- Almost colourless liquid with a characteristic fruity smell that becomes a vapour or gas once it has been poured into the Penthrox Inhaler

Terms & Conditions

Please note: In accordance with the NSW Ministry of Health's regulations, no in-office pickups are permitted. Penthrox products must be sent via courier directly to the licensed site.

We only supply to doctor's, dentists, paramedics and NSW first aiders with proper credentials. Your documentation will be verified before dispatch in accordance with the Poisons and Therapeutic Goods Act of 1966, as required by the Therapeutic Goods Administration of Australia.





PENTHROX INHALER – PACK 10

Brand: MDI

SKU: ME-759046-MK2

Description

Standard inhaler unit for single use. This product is designed to be used with the separate vial packs (MDI-ME-759044) as these units do not come with the methoxyflurane vial included.

Penthrox (methoxyflurane) is a non-narcotic inhalable analgesic that is self-administered by a conscious haemodynamically stable patient (under supervision of personnel trained in its use) via the hand-held Pentrox Inhaler to provide pain relief.

Standard Contents

- 10 x Pentrox Inhaler



PENTHROX 3ML VIAL – PACK 10

Brand: MDI

SKU: MDI-ME-759044

Description

Penthrox (methoxyflurane) is a non-narcotic inhalable analgesic that is self-administered by a conscious, haemodynamically stable patient (under supervision of personnel trained in its use) via the hand-held Pentrox Inhaler to provide pain relief.

Features

- It is clear
- Almost colourless liquid with a characteristic fruity smell that becomes a vapour or gas once it has been poured into the Pentrox Inhaler

Standard Contents

- 10 x 3ml Pentrox (methoxyflurane)
- Does not include inhalers (green whistle)



Penthrox training Kit with AC chamber

SKU: MDI-MEMS251



**MEGA
MEDICAL**

www.megamedical.com.au

Medical Equipment & Gases Australia

ABN 76 053 999 205

9 Middleton Road Cromer NSW 2099

1300 881 824 info@megamedical.com.au