



**MEGA
MEDICAL**

www.megamedical.com.au

**EMERGENCY OXYGEN &
RESUSCITATION + MEDICAL OXYGEN
CONSUMABLES + DEFIBRILLATORS**

Emergency Oxygen & Resuscitation	1
Medical Oxygen Masks	1
Nebulisers	2
Oxygen and Resuscitation Kits	3
Resuscitation Accessories	5
Medical Oxygen Consumables	8
Medical Oxygen Masks	8
Medical Oxygen Tubing	10
Accessories and Parts	10
Guedel Airways	10
Defibrillators	12

Introduction

Mega Medical is Australia's leading supplier of medical gas equipment, with more than 20 years of experience in the industry.

Australian-owned and operated, Mega Medical proudly serves all Australian states and territories health services and supplies medical equipment nationwide. We provide our customers with access to high-quality products, backed by the excellence of our customer sales and support division.



**MEGA
MEDICAL**

www.megamedical.com.au

EMERGENCY OXYGEN & RESUSCITATION

MEDICAL OXYGEN MASKS



RESUSCITATION MASK SILICONE - ADULT

Brand: Liberty

SKU: LR5134

Description

These masks are produced from medical-grade silicone, the patient mask comes in a paediatric size, suitable for children who may need more coverage. Offering a snug and comfortable fit, the mask helps teams to ensure effective treatment by providing a mask that suits the specific needs of each patient. Selecting the right sized mask helps to keep patients calm and relaxed during treatment, facilitating the best possible results from their care.

Variations



Resuscitation Mask
Silicone - Child
SKU: LR5133



Resuscitation Mask Silicone
- Large Adult
SKU: LR5135



RESUSCITATION MASK SILICONE – ADULT – (GALEMED)

Brand: Galemed

SKU: GA5134

Description

The silicone medical mask in adult size is designed for use with medical-grade oxygen and fits snugly and comfortably over the lower part of the patient's face.

Specifications

- Polycarbonate and/or silicone rubber
- Autoclavable up to 121°C
- Connections: 22mm

Variations



Resuscitation Mask
Silicone – Child –
(GaleMed)
SKU: GA5133



Resuscitation Mask
Silicone – Large
Adult – (GaleMed)
SKU: GA5135



NEBULISER COMPRESSOR

Brand: BestNeb

SKU: LRPNB

Description

Manufactured by BestNeb, this nebuliser compressor makes it easy for aerosol therapy to be delivered in a hospital, healthcare facility or home environment. The air and oxygen compressor is suitable for administration of medicated aerosol inhalation, something clinicians and patients may prefer to traditional medication delivery.

- Ultra quiet operation
- High quality aerosol compressor
- Large air intake filter
- Compartment for storing accessories

Comes complete with:

- 1 adult mask
- 1 child mask
- 1 nebuliser bowl
- 1 x 2-Metre supply of tubing
- 1 mouthpiece
- 1 T-piece
- 1 reservoir tube 15cm

Additional Purchase items:

- Spare filters – Pack of 5 (LRPNBF)

Features

- Material: Polycarbonate and/or silicone rubber
- Autoclavable up to 134°C
- Adjustable handles allows for easy movement.

Specifications

- Power: 220/240 V ~ 50 HZ
- Dimensions: 29 x 20 x 12 CM
- Weight: 2.5 Kg
- Warranty: 2 Years
- Maximum Pressure: 32 PSI
- Flow: 12 LPM
- Sound Level: 50 DB
- Particle Size (Mmad): 1.1 Micron



Child T-Piece Nebuliser Mask with 2.13m Tubing
SKU: BOS8906



Adult T-Piece Nebuliser Mask with 2.13m Tubing
SKU: BOS8924



Nebuliser Spare Filter – Pack of 5
SKU: LRPNBF



SMALL RESUSCITATION KIT – LITTLE RED

Brand: Mega Medical

SKU: mega002

Description

The 'Little Red' Medical Oxygen Resuscitation Pack is designed as a first response first aid kit for professionally trained medical personnel in the pre-hospital environment.

Meets Australian Standard AS2488 – 1995
This product is entered in the ARTG 155544

Features

- Easy to clean, durable lightweight material plus highly reflective tape for safety
- Comfortable handles and shoulder strap for ease of carrying
- Sleeve for securing the medical oxygen bottle
- Maintenance repairs and preventative maintenance programs to meet AS. 3551 – Available from Mega.

Standard Contents

- 'Little Red' Resuscitation Bag
- Medical Oxygen Regulator/Flowmeter 1-25LPM
- Disposable Silicone Adult Resuscitator
- Re-Useable Silicone Child Resuscitator Mask
- Bacterial Filter
- Adult Oxygen Mask with Tubing
- Child Oxygen Mask with Tubing
- Airways 100mm (Red)
- Airways 90mm (Yellow)
- Airways 80mm (Green)
- Medical Oxygen Resuscitator Instruction Manual



Compatible Cylinder Sizes

- C and CH Size standard Cylinder
- C – IVR Cylinder (Coregas)
- C – Inhalo Cylinder (BOC Gases)
- Take O2 2.0 (Air Liquide)

Related Products



Small Resuscitation Kit – Little Red – Without Regulator
SKU: mega003



LARGE RESUSCITATION KIT – BIG RED



Brand: Mega Medical

SKU: mega004

Description

The 'Big Red' medical oxygen resuscitation back pack with a defibrillator pouch is designed as a first response first aid kit for the professionally trained medical person in the pre-hospital environment. Having a durable, light-weight concealed backpack harness allows for hands-free movement in open and confined spaces. A variety of different sized internal pouches gives many options for accessories including optional suction and/or demand valves.

Meets Australian Standard AS2488 – 1995
This product is entered in the ARTG 155544

Specifications

- Weight (without Cylinder): Approx 2.5kg
- Weight (with C size Cylinder): Approx 6.5kg
- Length: 620mm
- Width: 320mm
- Diameter: 120mm
- Cylinder Holder Size: C Size – 410 Litres

Features

- Easy to clean durable lightweight material plus highly reflective tape for safety
- Comfortable side handles for ease of carrying
- Defibrillator pouch fits all modern defibrillators
- Variety of internal pocket sizes with clear windows
- Concealed backpack straps/harness
- Reinforced base and sleeve for securing the medical oxygen bottle
- Maintenance repairs and preventative maintenance programs to meet AS. 3551 – Available from Mega.

Related Products



Large Resuscitation Kit – Big Red – With Suction
SKU: mega005



Combat Sports Resuscitation Kit
SKU: Mega-Combat-Sports-Kit-Ex



Emergency Oxygen Kit – Equipment Only
SKU: mega-emerg-no-cylinder



OXYGEN DEMAND VALVE RESUSCITATOR + 1M HOSE AND MASK

Brand: Mega Medical

SKU: MEGA-DVR

Description

The demand valve resuscitator provides a safe and effective means of delivering an oxygen concentration to ventilate non-breathing or respiratory distressed patients. The unit can also be used as a standard demand valve to deliver oxygen to spontaneously breathing patients upon demand. This is quality medical oxygen equipment built by Ambulac with an in-built oxygen flow control valve.

Features

- Safe and effective means of delivering oxygen
- Can be used as a demand valve resuscitator
- Integral pressure relief valve
- Can be attached to any standard 400 kPa oxygen outlet
- Includes 1m hose assembly and LR 5134 reusable mask silicone

Specifications

- Manual Flowrate (LPM): 40 – 60
- Demand Flowrate (LPM): 0 – 160 LPM
- Demand Flow Triggering System: -5cm h30 @ 100 LPM, maximum
- Supply Pressure (kPa): 400 nominal
- Pressure Relief Valve: 54 – 60 cm h30
- Exhalation Resistance: -5 cm h30 @ 50 LPM, maximum
- Inlet Connection: 9/16" – 18 D.I.S.S
- Patient Connection: 22 / 15mm male / female taper
- Weight (grams): 230

Related Products



Demand Resuscitator – Diaphragm Valve
SKU: CIG552082



Patient Valve Assembly
SKU: CIG552083



Oxygen Demand Valve Resuscitator
SKU: AMB-DVR



Oxygen Hose Assembly – Demand Fitting – 1m
SKU: CIGMA169



Pocket Resuscitation Mask with Filter
SKU: LR1325



Fingertip Pulse Oximeter – Oxygen Saturation Monitor
SKU: FPO100B



ADULT RESUSCITATOR WITH POP-OFF – DISPOSABLE

Brand: Besmed

SKU: 55786

Description

The adult resuscitator is for single use (non autoclavable) and is made from high quality PVC and silicone. The unit is suitable for use as a manual resuscitator on adult patients and is made by Besmed which is one of the best known brands on the market.

Specifications

- Material: Polycarbonate, silicone rubber and polyvinyl chloride
- Connectors: 15/22mm ISO patient connector, 23mm ID silicone neck, 25mm ID (to silicone bag inlet), 25mm OD (Intake valve port), 6mm OD (supp gas inlet)
- Inspiratory / Expiratory Resistance: <5 cm h30 at 50 LPM flow
- Reservoir Volume: 2500 mL
- Max. Oxygen Concentration: 99% (with oxygen supply)
- Operating Condition: -18°C to 50°C
- Stroke Volume: One Hand = 700ml Two Hands = 900ml

Features

- Unit is perfect for use in emergency oxygen kits
- Packs away for efficient storage

Related Products



Infant Resuscitator – Disposable
SKU: LR2297



Resuscitator Disposable – Child
SKU: LR2287



Silicone Adult Resuscitator – Reusable
SKU: LRES3103RI



Silicone Child Resuscitator Reusable
SKU: LR1406



Oxygen Reservoir Adult 2500 ML – Single Use
SKU: LR24005



Bacterial Filter
SKU: LR70015



HANDHELD MANUAL SUCTION PUMP

Brand: Mega Medical

SKU: NAN-4217

Description

The Handheld Manual Suction Pump features a simple and durable design that makes it easy to operate with one hand. An adjustable pressure knob provides 2 settings – 50% or 100% (100% provides a maximum vacuum suction of 39.9 kPa/300mmHg), this reusable device is perfect for suctioning on the go or any other situation where power may not be available. Its ergonomic handle provides a comfortable and secure grip, while the compact size makes it easy to store and transport.

Specifications

- Vacuum at 50% – 20kPa
- Vacuum at 100% – 39.9kPa
- Canister volume – 300mL
- Total weight – 230g

The Handheld Manual Suction Pump is made of high-grade plastic which can easily be cleaned between uses. The pump comes with 2 different catheter sizes for pharyngeal suction of adult and paediatric patients in emergency situations.

Medical Oxygen Consumables

MEDICAL OXYGEN MASKS

ADULT NASAL CANNULA WITH 2M TUBING



Brand: Salter Labs

SKU: LR3233

Description

The nasal cannula is used to deliver supplemental oxygen to a patient in need of respiratory help. Also referred to as 'Nasal Prongs' these are used to administer oxygen to the patient via the nose to increase efficiency in gas delivery and provide a more comfortable and convenient delivery method.

Specifications

- Length: 2m
- Colour: Clear
- Material: Medical PVC

Features

- Soft, transparent, odour free vinyl to ensure the best fit and user comfort.
- Flexible oxygen tubing prevents occlusion when tubing is kinked.

Standard Contents

- Soft tip adult nasal cannula
- Patient cautionary notes

Patient Notes: This product is designed for single patient use and is not intended to be reprocessed. Reusing, soaking, rinsing or sterilising this device with chemicals may leave harmful residues and may impede the functions of this device. Reusing single patient use device may also increase the risk of cross contamination. It is strongly advised that the product be used according to trained medical personnel's guidance and the instructions for use before the marked expiration date.

Related Products

50

Adult Nasal Cannula with 2m Tubing - Box 50
SKU: LR3233-Box50



Child Nasal Cannula with 2m tubing
SKU: LR3232

50

Child Nasal Cannula with 2m tubing - Box 50
SKU: LR3232-Box50



OXYGEN MASK WITH 2M TUBING – ADULT

Brand: Salter Labs

SKU: LR3645

Description


Adult oxygen therapy mask or 'Hudson Mask' as it is commonly referred to, is the ideal solution in patient oxygen therapy. This oxygen therapy mask is manufactured by Salter Labs and is one of the premium masks on the market. Moulded from superior odour-free vinyl, this mask is made for ultimate function and comfort. The mask features a thick elastic strap and an adjustable aluminium nose clip to ensure the most secure and comfortable fit for the patient. The 2 metre translucent tubing will connect to all medical oxygen therapy units with an oxygen outlet nipple.

Features

- Soft, transparent, odour-free vinyl
- Aluminium fixer and elastic strap to ensure the best fit and user comfort

Related Products

50



Oxygen Mask with 2m Tubing - Adult - Box 50
SKU: LR3645 - Box50



Oxygen Mask with 2m Tubing – Child
SKU: LR3644

50



Oxygen Mask with 2m Tubing - Child - Box 50
SKU: LR3644-Box50



Child High-Oxygen Concentration Mask
SKU: LR3674



Adult Oxygen Concentration Non-Rebreathing Mask
SKU: LR3675

MEDICAL OXYGEN TUBING



Oxygen Tubing 10m
SKU: LR31610



Oxygen Tubing – 2m
SKU: LR2002



Oxygen Tubing Connector
SKU: BOS1215

ACCESSORIES AND PARTS



Oxygen Humidifier Bottle
Disposable
SKU: BOS7600



Oxygen Tubing Y Connector /
Flow Splitter
SKU: PAAPM88409

GUEDEL AIRWAYS



Standard Contents

- One (1) of each
- 60mm – Black
 - 70mm – White
 - 80mm – Green
 - 90mm – Yellow
 - 100mm – Red
 - 110mm – Light Blue

GUEDEL AIRWAYS – SET OF 6

Brand: Mega Medical

SKU: LR6300

Description

The Guedel Airway is designed for free air flow by protecting from airway occlusion from the relaxation of the tongue in an unconscious patient. The Guedel Airway devices have a central lumen which allows for various suction catheters to be inserted.

The top section is colour coded for rapid identification of size as well as a thickened area for protection from biting. This bite block is ideal with seizures and unprotected airways. It also has benefits for the intubated patient preventing them from crushing the endotracheal tube.

Features

- Consists of a large flange, a reinforced bite portion and a gentle curve that follows the contour of the tongue
- Colour code for size identification conforms to ISO requirements
- Meets ISO-10993:1997 for biocompatibility requirements



GUEDEL AIRWAY – 40MM – PINK

Brand: Mega Medical

SKU: LR6304

Description

The Guedel Airway is designed for free air flow by protecting from airway occlusion from the relaxation of the tongue in an unconscious patient. The Guedel Airway devices have a central lumen which allows for various suction catheters to be inserted.

The top section is colour coded for rapid identification of size as well as a thickened area for protection from biting. This bite block is ideal with seizures and unprotected airways. It also has benefits for the intubated patient preventing them from crushing the endotracheal tube.

Features

- Consists of a large flange, a reinforced bite portion, and a gentle curve that follows the contour of the tongue
- Colour code for size identification conforms to ISO requirements
- Meets ISO-10993:1997 for biocompatibility requirements



Variations



Guedel Airway
50mm – Blue
SKU: LR6305



Guedel Airway
80mm – Green
SKU: LR6308



Guedel Airway
100mm – Red
SKU: LR6310



Guedel Airway
70mm – White
SKU: LR6307



Guedel Airway
60mm – Black
SKU: LR6306



Guedel Airway
110mm – Light Blue
SKU: LR6311



Guedel Airway
90mm – Yellow
SKU: LR6309

Boxes Of 50



Guedel Airway - 60mm
- Black - Box 50
SKU: LR6306-Box 50



Guedel Airway - 90mm
- Yellow - Box 50
SKU: LR6309-Box 50



Guedel Airway - 80mm
Green - Box 50
SKU: LR6308-Box 50



Guedel Airway - 70mm
White - Box 50
SKU: LR6307-Box 50



Guedel Airway -
100mm - Red - Box 50
SKU: LR6310 - Box 50



Guedel Airway
110mm - Light Blue
Box 50
SKU: LR6311 - Box 50

Defibrillators



HEARTSINE SAMARITAN PAD 360P FULLY AUTOMATIC AED

Brand: Heartsine

SKU: PAD-360P

Description

The HeartSine 360P Defibrillator (AED) is a device designed to save lives with its fully automatic shock delivery system. With its user-friendly interface, even non-medical personnel can confidently operate it. Its safety features ensure 100% effectiveness in delivering visual and verbal prompts to guide rescuers through CPR procedures. The unit weighs only 2kg making it lightweight and portable, and its high Ingress Protection rating of IP56 allows it to function in outdoor and harsh weather conditions.

Features

- 8-year warranty
- Automatic Shock Delivery
- Easy to Follow Visual and Verbal Guides
- Durable
- Advanced Technology
- Single Battery/ Electrode Cartridge

Includes:

- HeartSine Samaritan PAD 360P Defibrillator
- Protective cover
- Defibrillator Pads
- Defibrillator Battery
- 8 Year Warranty

Related Products



Heartsine Paediatric Pad-Pak – Pads & Battery
SKU: PAD-PAK-04



Heartsine Adult Pad-Pak – Pads & Battery
SKU: PAD-PAK-03



Heartsine Samaritan PAD 500P Semi-Automatic AED-
SKU: PAD-500P



Heartsine Samaritan PAD 350P Semi-Automatic AED-
SKU: PAD-350P



LIFEPAK CR2 ESSENTIAL FULLY-AUTOMATIC AED

Brand: LifePak

SKU: CR-2E-A

Description

The LIFEPAK CR2 Essential Fully Automatic defibrillator is the ideal solution for those seeking guaranteed performance during critical moments. With the powerful AED Management Software – LIFELINKcentral – this defibrillator is always ready to save lives.

Monitor with Ease

Thanks to LIFELINKcentral, users can monitor every CR2 connected to the software via Wi-Fi. This software generates reports and alerts users automatically if any connected AEDs are not ready, such as when the battery is not charged.

Features

- 4 years of standby battery life
- Durable
- 8-year warranty
- Weight: 2kg
- Dimensions: 34 × 31 × 20 cm

Related Products



LifePak CR2 Battery
SKU: CR-2-BAT



LifePak CR2 Replacement Electrode Kit
SKU: CR-2-PAD



LifePak CR2 Fully Automatic AED – WiFi
SKU: CR-2-A



LifePak CR2 Essential Semi-Automatic AED
SKU: CR-2E



LifePak CR2 Semi-Automatic AED – WiFi
SKU: CR-2



AED Basic Prep Kit
SKU: CA-PK2



AED Premium Prep Kit
SKU: CA-PK1



CardiAct Poly AED Wall Sign-
SKU: CA-07



CardiAct Standard AED Cabinet
SKU: CC-25



**MEGA
MEDICAL**

www.megamedical.com.au

Medical Equipment & Gases Australia

ABN 76 053 999 205

9 Middleton Road Cromer NSW 2099

1300 881 824 info@megamedical.com.au